



HAORUN MEDICAL PRODUCTS CO.,LTD



CONTENTS

page 01~02	Oral/Enteral Syringes Series
page 03~04	Disposable Syringes Series
page 05	Infusion Sets Series
page 06	Blood Transfer Series
page 06	Medical Staplers Series
page 07	Disposable Needles Series
page 08	Masks Series
page 09	Urine/Stool Containers Series
page 10	Medical Tapes Series
page 10	Medical Gloves Series
page 10	Hernia Mesh Series

Oral/Enteral Syringes Series



Oral/Enteral Syringe with Tip Cap

Item code: CMOSC

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- **Non-luer tip, can not fit with hypodermic needles**
- Latex free, medical grade PP
- Transparent barrel, colored plunger
- **With tip cap for racap after use**



Oral/Enteral Syringe with Adaptor

Item code: CMOSA

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- **Non-luer tip, can not fit with hypodermic needles**
- Latex free, medical grade PP
- Transparent barrel, colored plunger
- **With transparent adaptor for bottle, fitting the bottle neck ID 13mm~23mm**



Oral Syringe with Dispensing Tube

Item code: CMODT

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- **Non-luer tip, can not fit with hypodermic needles**
- Latex free, medical grade PP
- Transparent barrel, colored plunger
- **With dispensing tube (color and length available)**



Oral/Enteral Syringes in Packages



**Oral Syringe (Luer Slip Tip)
with Dispensing Tube**

- Item code: CMOSDT(LS)*
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
 - Luer tip (Hypodermic needles be available)
 - With protective cap and dispensing tube



**Oral Syringe (Two Part)
with Adaptor**

- Item code: CMOSA(2P)*
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
 - Luer tip (Hypodermic needles be available)
 - With adaptor



Irrigation Syringes

- Item code: CMSC*
- Size: 20ml, 30ml, 50ml, 60ml, 100ml, 150ml
 - 50ml with bulb



**Oral/Enteral Syringe Female
Reverse Luer Lock**

- Item code: CMOSLL*

Disposable Syringes Series



Disposable Three-parts Syringe Luer Slip

Item code: CMSS

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml
- Luer slip, with or without needle



Disposable Three-parts Syringe Luer Lock

Item code: CMSL

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml, 100ml
- Luer lock, with or without needle



Disposable Two-parts Syringes

Item code: CMST

- Size: 2ml, 3ml, 5ml, 10ml, 20ml
- Luer slip, with or without needle



Auto Disable Syringe

Item code: CMSAD

- Size: 0.05ml, 0.5ml, 1ml, 2/3ml, 5ml, 10ml
- Simple hypodermic needle available
- **Plunger be locked in the barrel after use**



Safety Syringe with Retractable Needle (Manually)

Item code: CMSA

- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Simple hypodermic needle available
- **Needle be retractable into the barrel and locked after use**



Insulin Syringe

Item code: CMSI


- Size: 0.3ml, 0.5ml, 1ml
- Fixed needle, 27g, 28g, 29g, 30g, 31g

Infusion Sets Series



Precision Infusion Set

Item code: CMIF-JM



Lightproof Infusion Set

Item code: CMIF-LP



Infusion Set

Item code: CMIF-BR



Burette Infusion Set

Item code: CMIFB

• Size: 100ml, 150ml • 60 drops = 1 ± 0.1 ml



Item code: CMIF-01



Item code: CMIF-02



Item code: CMIF-03



Item code: CMIF-04

Disposable Infusion Set / I.V. Set.

Material: medical grade PVC, with or without DEHP

Specification: 1.big chamber with air vent 2.common pvc or elastic tube, 120~180mm
3.rubber or Y injection site 4.polybag or blister packing



Scalp Vein Set

Item code: CMIV

• luer slip or luer lock • size: 16G~27G



I.V. Cannula

Item code: CMCI-2

• Size: 14G, 16G, 17G, 18G, 20G, 21G, 22G, 24G

Blood Transfer Series



Blood Bag (Triple)

CPD-63ml for collection of 450ml blood & SAG-M solution 100ml & transfer bag



Blood Bag (Double)

CPDA-1-63ml for collection of 450ml blood & transfer bag



Blood Bag (Single)

CPDA-1-63ml for collection of 450ml blood



Blood Transfusion Set

Item code: CMBT

Medical Staplers Series



Disposable Linear Stapler

Item code: CMMS-01



Disposable Linear Cutter Stapler

Item code: CMMS-02



Disposable Circular Stapler


Item code: CMMS-03



Disposable Hemorrhoids Stapler

Item code: CMMS-04


Disposable Needles Series



Hypodermic Needle

Item code: CMHN

Size: 14G, 15G, 16G, 17G, 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G




Spinal Needle

Quincke Bevel & Pencil Point

Gauge: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G

Specification: 90mm



Special Needles

Customers From All Over The World





HAORUN MEDICAL PRODUCTS CO.,LTD

TEL 0086-574-87706193

WHATSAPP/WECHAT/TEL:008613566542343

EMAIL:SALES@HAORUNMED.COM

WWW.HAORUNMED.COM