



**HAORUN MEDICAL PRODUCTS CO.LTD**



# Safety Double Lumen Laryngeal Mask Silicone

## Characterics



- Safer structure**  
 New structure with double reflux drainage channels to prevent reflux aspiration and prevent complications effectively
- Multi-function**  
 Bronchofibroscope and endotracheal tube could be used together
- Good ventilation performance**  
 Good sealing and ventilation effect
- High quality material**  
 Made by medical elastic material, soft and sleek, less damage to oral soft tissue
- More comforts**  
 Special design to avoid nerve ending areas, more comfortable for patients.

## Double reflux fluid drainage channels

Reflux fluid drainage channel one



Reflux fluid drainage channel two



Suction pipe to drainage reflux fluid

## Models & Package

Model	Sizes	Weight	Maximum Inflation	Packing
DM- H10	1.0#	Infant < 5 kg	4 ml	5 or 10 pcs/box , 50 pcs/ carton
DM- H15	1.5#	Infant 5- 12kg	7 ml	
DM- H20	2.0#	Child 12-25kg	10 ml	
DM- H25	2.5#	Child 25-40kg	14 ml	
DM- H30	3.0#	Adult 40- 65kg	20 ml	
DM- H40	4.0#	Adult 65- 90kg	30 ml	
DM- H50	5.0#	Adult > 90 kg	40 ml	

# Safety Oropharyngeal Airway



## • More comforts

- Design avoids the uvula and the central region of the tongue, less stimulation
- With lip protection sheet
- Made by medical elastic material with super slippery surface, soft and comfortable

Size: Pediatric / Adult Normal / Adult Big




## • Special design, safer for patients


- With independent endotracheal tube channel & independent ventilation channel
- Can hold the tongue, prevents airway obstruction by tongue falls back.

## • the anesthesia operation process simplified

### • Current anesthesia process



### • Anesthesia process by using DM safety oropharyngeal airway



*Apply to whole anesthesia process*

## Oropharyngeal Guedel Airway



- Made by medical PVC
- Atraumatic rounded tip
- With reinforcement insert
- Color coding for size identification
- Latex free, sterile by EO, single use
- Individual packaging




Available size: 0# 1# 2# 3# 4# 5# 6# 7# 8# 9# 10# 11# 12#

# Bite Block



## Lip protective type

Size: Pediatric / Adult Normal / Adult Big



## Smart type



Size: Pediatric / Adult Normal / Adult Big

- **Tube protection channel**  
Fix and protect the endotracheal tube
- **Tube fixed cylinder:**  
The endotracheal tube can be wound on the cylinder to be fixed
- **Suction channel**  
Assist sputum suction
- **Medical elastic material:**  
Soft, more comfortable, better for teeth

### *Multi-function bite block DAMI Brand*

Endotracheal tube protection function + Sputum aspiration function + Tooth protection function

# One-hand-hold Anesthesia Mask



## Characterics



- **Innovative craft air cushion**

Innovative blow molding process production brings the soft air cushion, which fits the face better and seals well, avoids air leakage.

- **Innovative one-hand-operation structure**

The mask can be grasped with one hand and the finger fatigue caused by long-term holding can also be avoided.

Size: 3# 4# 5# 6#



New developed



“Finger buckle” hold method — Able to be hold with one hand easily

---

# High Quality Air Cushion Anesthesia Mask



## Characterics

- Anesthesia mask is perfectly designed according to facial structure
- The seal is soft to provide maximum comfort
- non-toxic, harmless, high strength
- Wave shape buffer

Size: 1# 2# 3# 4# 5# 6#

# Safety Oxygen Mask

New coming: Will be on the market soon

A new generation oxygen mask



**Multi-function oxygen mask**  
**All round protect for medical workers**  
**More comfortable for patients**

There is an opening with a flap design in the middle can be used for painless diagnosis and treatment, intubation and other operations

Both sides are designed with carbon dioxide outlet and filter element, which can effectively prevent droplet and aerosol spreading

**Size:** Pediatric / Adult

### Application:

- ① Anesthesia recovery room: For oxygen supply & Protect patients and medical personnels.
- ② Endoscopy center
- ③ Infection, ICU, respiratory and other departments

# PVC Oxygen Mask

- . The star lumen tubing can ensure oxygen flow even if the tube is kinked, different length of the tubing is available
- . Material: Clear, soft PVC for patient comfort and visual assessment.
- . The product can be transparent green and transparent white.



Adjustable nose clip assures comfortable fit.

The turn up rim ensure comfortable fit with good seal.



# Endotracheal Tube



Endotracheal Tube Normal with Cuff	
Endotracheal Tube Reinforced with Cuff	
ID 2.0mm	ID 6.0
ID 2.5mm	ID 6.5
ID 3.0	ID 7.0
ID 3.5	ID 7.5
ID 4.0	ID 8.0
ID 4.5	ID 8.5
ID 5.0	ID 9.0
ID 5.5	ID 9.5

## Guided Suction Catheter

New coming: Will be on the market soon



- **Shaping function**

After inserting the suction catheter into the endotracheal tube, the endotracheal tube can be easily shaped according to the position of the patient's glottis.

- **Designed with an individual breathing channel**

It can connect with oxygen supply device





## **HAORUN MEDICAL PRODUCTS CO.LTD**

- Safer
- More Comfortable
- More convenient to

Tell: 0086-574-87706193

Web: [www.haorunmed.com](http://www.haorunmed.com)